

## **Center Scripts (center\_scripts\_v1.0) Walk-Through**

### **July 31, 2007**

#### **Introduction**

Center scripts are used to do SPM batch analysis. This version of center scripts is compatible with SPM5 version. The following things must be done before running the analysis:

- Add center\_scripts\_v1.0 directory and its sub-directories on MATLAB path.
- Data must be organized such that each subject has run directories.
- All the required parameters for the analysis must be entered in a preference file. A detailed explanation of the preference file is given in the next section. Example preference files like 'cs\_prefs\_aod\_analyze.m', 'cs\_prefs\_aod\_nifti.m' and 'cs\_prefs\_structural.m' are located in the directory center\_scripts\_v1.0/prefs.
- Create your own 'spm\_defaults.m' file. Example 'spm\_defaults.m' file is located in the directory center\_scripts\_v1.0.

#### **Note:**

1. Preference and spm\_defaults files must be added to the MATLAB path.
2. Use 'spm\_defaults.m' file in directory center\_scripts\_v1.0 instead of the standard 'spm\_defaults.m' file as it contains additional fields that are required for proper running of the center scripts.

You can run the analysis using two functions like 'cs\_run\_all' and 'cs\_run\_sub'. The syntax for each function is given below:

- 'cs\_run\_all' – This function is used for analyzing specified subject directories. The syntax for running this file is as follows:
  - `prefs_file='cs_prefs_aod_nifti';subject_directory='J:\AOD_Data\2subject s-unprocessed\s01\';`
  - `cs_run_all(prefs_file, subject_directory);`
- 'cs\_run\_sub' – This function is used to automatically get the subject directories based on field scandir\_regexp in structure csprefs (Global variable in preference file). It can also be used if you have several tasks to be run at once. Place each task in a preference file and call at the MATLAB command prompt as follows:
  - `prefs_files = {'cs_prefs_aod_nifti', 'cs_prefs_structural'};`
  - `cs_run_sub(prefs_files);`

## Preferences file

Preferences file contains the parameters required to do pre-processing or stats. To run a particular step use a value of 1 and a value of 0 means don't run that step. The following are the steps involved:

## Processing Steps

The following are the main steps involved:

- a. Dicom Converter (csprefs.run\_dicom\_convert)
- b. Realign (csprefs.run\_realign)
- c. Co-register (csprefs.run\_coregister)
- d. Slice Time correction (csprefs.run\_slicetime)
- e. Normalize (csprefs.run\_normalize)
- f. Smooth (csprefs.run\_smooth)
- g. Filter (csprefs.run\_filter)
- h. Stats (csprefs.run\_stats)
- i. Segment (csprefs.run\_segment)

## Data directory

Data directory and log file information is described below:

- a. csprefs.exp\_dir – Experimental directory or the root directory where all the data is located.
- b. csprefs.logfile – Log file containing information whether center scripts ran successfully or not.
- c. csprefs.errorlog – Error information is stored in this log file.
- d. csprefs.spm\_defaults\_dir – Enter full file path of the copy of 'spm\_defaults.m' file.
- e. csprefs.scandir\_regexp – Regular expression for subject directories. This field will be used by the 'cs\_run\_sub' function to get the subject directories. Some of the standard regular expressions are described in the Regular Expressions section.
- f. csprefs.rundir\_regexp - Regular expression for run directories.
- g. csprefs.scandir\_postpend - Regular expression used to add an additional path to a subject directory. For example if 'C:\sub01\' is the subject directory and you have run directories inside 'Study001' directory, regular expression will be 'Study\d{3}'. Leave csprefs.scandir\_postpend as empty if you have run directories immediately inside a subject directory.
- h. csprefs.rundir\_postpend - Regular expression used to add an additional path to a run directory. For example if 'C:\sub01\Study001\aod\_v1\_001' is the run directory and you have images inside 'Original\Nifti' directory, regular expression will be 'Original\\Nifti' ('Original\Nifti' under Unix). Leave csprefs.rundir\_postpend as empty if you have images immediately inside a run directory.
- i. csprefs.dummyscans – No. of dummy scans. All the dummy scans will be moved to directory 'dummies'. Leave this as 0 if you want to use all the scans. In case of

- a 4D Nifti file, all the dummy scans will be written with the same file name in 'dummies' directory.
- j. csprefs.tr – TR for your experiment.

## **Dicom Converter**

Variables contained in dicom converter as follows:

- a. csprefs.dicom.file\_pattern – File pattern of dicom files.
- b. csprefs.dicom.format - Options for converting dicom files.
  - a. '3d\_analyze' – Dicom files will be converted to 3D analyze images.
  - b. '3d\_nifti' – Dicom files will be converted to 3D Nifti images.
  - c. '4d\_nifti' – Dicom files will be converted to a 4D Nifti image.
- c. csprefs.dicom.write\_file\_prefix – Prefix for the image files that are converted from dicom files.
- d. csprefs.dicom.outputDir – Output directory to write analyze or Nifti images.

## **Realign**

Variables contained in realign are as follows:

- a. csprefs.coregister – Coregister images to the subject's first scan first session. If you have already performed this step, you need to set this field as 0.
- b. csprefs.reslice – Reslice images using 'spm\_reslice' function.
- c. csprefs.use\_inrialign – Uses INRIAlign function to realign. A value of 0 means 'spm\_realign' function is used.
- d. csprefs.realign\_pattern – File pattern of the images that need to be realigned.
- e. INRIAlign defaults:
  - a. csprefs.inrialign\_rho – Rho function for INRIAlign function. Default value is 'geman' (Geman-McClure function). Other options are as follows:
    - i. 'absolute' - Quite slow and not very robust
    - ii. 'huber' - Huber function
    - iii. 'cauchy' - Cauchy function
    - iv. 'leclerc' - Leclerc-Welsch function
    - v. 'tukey' - Tukey's biweight function
  - b. csprefs.inrialign\_cutoff – Cut off distance. Default is 2.5.
  - c. csprefs.inrialign\_quality – Quality. Default is 1 (slow and high quality).
- f. csprefs.realign\_fwhm – Size of smoothing kernel. Default is 8.
- g. csprefs.realign\_rtm – Realign to the mean image. Default is 0. A value of 1 does not work for INRIAlign.
- h. csprefs.reslice\_write\_imgs – Option for writing re-sliced images. Default is 0.
- i. csprefs.reslice\_write\_mean – Option for writing mean image. Default is 1.

## **Co-register**

Variables contained in co-register step are as follows:

- a. `csprefs.run_coreg` – Run co-register step. Options are 1 and 0. A value of 1 will estimate the registration parameters using the reference image (`csprefs.coreg.ref`) and the source image (`csprefs.coreg.source`).
- b. `csprefs.run_reslice` – Run re-slice step. Options are 1 and 0. 1 means images will be re-sliced using the reference image (`csprefs.write.ref`). The images that will be re-sliced are source image (`csprefs.coreg.source`) and the other images (images obtained using file pattern from variable `csprefs.coreg.other_pattern`). The new set of images will have prefix 'r'.
- c. `csprefs.coreg.ref` – Reference or template image used for registration.
- d. `csprefs.coreg.source` – Source image or file pattern.
- e. `csprefs.coreg.other_pattern` – You can specify additional images to apply the registration parameters.
- f. `csprefs.coreg.write.ref` – Reference image used for re-slicing source and other images.

**Note:** After the co-register step, the headers of the images will be modified to incorporate the registration parameters.

### **Slice Time Correction**

Variables contained in slice time are as follows:

- a. `csprefs.slicetime_pattern` – File pattern of the images that need to be slice time corrected.
- b. `csprefs.sliceorder` – Slice time order. Include all the slices.
- c. `csprefs.refslice` – Reference slice. Default is middle slice.
- d. `csprefs.ta` – Time of acquisition. Leave for now this variable on 'default'.

### **Normalize**

Variables contained in normalize are as follows:

- a. `csprefs.determine_params` – Determine normalization parameters. If you have already determined normalization parameters, set this value to 0 and specify the MAT file in field `csprefs.writenorm_matname`.
- b. `csprefs.write_normalized` – Option for writing normalized images. Set this value to 1 for now.
- c. `csprefs.params_template` – Template image used for parameter estimation. Default is 'EPI.nii'.
- d. `csprefs.params_pattern` – Name of image to be used for parameter estimation. Usually this is the mean image created during realignment.
- e. `csprefs.writenorm_pattern` – File pattern of the images that require the parameters to be applied.
- f. `csprefs.writenorm_matname` - File or filter pattern containing the normalized parameters. To use this file you need to set field `csprefs.determine_params` to 0.

## Smooth

Variables contained in smoothing are as follows:

- a. `csprefs.smooth_kernel` – Smoothness kernel to be applied. Default is [10,10,10].
- b. `csprefs.smooth_pattern` – File pattern of the images that need to be smoothed.

## Filter

Variables contained in filtering are as follows:

- a. `csprefs.filter_pattern` – File pattern of the images that need to be filtered.
- b. `csprefs.cutoff_freq` – Normalized cut-off frequency. Default is 0.25.

## Stats

Variables in stats are as follows:

- a. `csprefs.stats_make_asciis` – For now leave this as 0.
- b. `csprefs.stats_ascii_script` – Full file path of the ascii script if variable `csprefs.stats_make_asciis` is set to 1.
- c. `csprefs.stats_beh_dir_name` – Behavioral directory name. This will be used only when make ascii script is used.
- d. `csprefs.stats_files_relative_path_sub` – Option is provided to find the onset and duration files relative to a subject directory or corresponding run directory. A value of 1 means onset files will be searched relative to a subject directory and a value of 0 means files will be searched relative to corresponding run directory.
- e. `csprefs.stats_dir_name` – Center scripts will create a directory using this name for storing the stats results.
- f. `csprefs.stats_pattern` – Stats will be done using this image file pattern.
- g. `sprefs.stats_beh_units` – Behavioral units. Options are 'scans' and 'secs'.
- h. `csprefs.stats_volterra` – Model interactions. For now leave it as 0.
- i. `csprefs.stats_basis_func` – Basis function to use. Options are 1 and 2.
  - a. 1 – 'hrf'
  - b. 2 – 'hrf (with time derivative)'
- j. `csprefs.stats_onset_files` – Specify onset files location in a cell array like `{'beh/TRG_PR_1_noslice.asc','beh/NOV_OM_1_noslice.asc','beh/STD_OM_1_noslice.asc';'beh/TRG_PR_2_noslice.asc','beh/NOV_OM_2_noslice.asc','beh/STD_OM_2_noslice.asc'}` where rows indicate onset files for that session. You can enter full file path if the onsets information doesn't vary between the subjects.
- k. `csprefs.stats_duration_files` – Number of files must match the number of onset files. For short events enter 0.
- l. `csprefs.stats_time_modulation` – To run time modulation you need to specify a cell array of size number of sessions by number of conditions. If you are not running time modulation, set variable `csprefs.stats_time_modulation` as 0 or empty cell array. In the example below we are using first order time modulation for two sessions and three conditions:

```
csprefs.stats_time_modulation = {1, 1, 1; 1, 1, 1};
```

- m. `csprefs.stats_parametric_modulation` - Parametric modulation can be specified for each condition each session. If you are not running parametric modulation, set variable `csprefs.stats_parametric_modulation` as 0 or empty cell array. In the example below we are using a first order polynomial for two sessions and three conditions:

```
csprefs.stats_parametric_modulation = {{'Targets', [1:23], 1}, {'Novels', [1:23], 1}, {'Standards', [1:184], 1}; {'Targets', [1:24], 1}, {'Novels', [1:23], 1}, {'Standards', [1:185], 1}};
```

- n. `csprefs.stats_global_fx` – Remove global effects. Options are 1 and 0. For now leave it as 0.
- o. `csprefs.stats_highpass_cutoff` - Number of seconds for high-pass filter. Default is 128. Put Inf (without quotes) for no filtering.
- p. `csprefs.stats_serial_corr` – Correct for serial corrections. Options are 1 and 0. For now leave it as 0.
- q. `csprefs.stats_tcontrasts` – Contrasts are specified in a matrix whose dimensions are number of contrasts by number of regressors. Number of regressors is equal to (number of sessions \* number of basis functions \* number of conditions \* (1 + order of time modulation + polynomial order of parametric modulation)) + number of sessions.
- r. `csprefs.stats_tcontrast_names` – Contrast names must be entered in a cell array whose length must equal the number of rows in variable `csprefs.stats_tcontrasts`.

## Segment

Variables in segmentation step are as follows:

- a. `csprefs.segment.pattern` – File pattern of the images that need to be segmented.
- b. The options for writing the output grey matter, white matter and CSF images are as follows:
- `[0 0 0]` - None
  - `[0 0 1]` - Native Space
  - `[0 1 0]` - Unmodulated Normalised
  - `[1 0 0]` - Modulated Normalised
  - `[0 1 1]` - Native + Unmodulated Normalised
  - `[1 0 1]` - Native + Modulated Normalised
  - `[1 1 1]` - Native + Modulated + Unmodulated
  - `[1 1 0]` - Modulated + Unmodulated Normalised
- c. `csprefs.segment.output.GM` – Option for writing the grey matter. Default value is `[1, 1, 1]`.
- d. `csprefs.segment.output.WM` – Option for writing the white matter. Default value is `[0, 0, 1]`.
- e. `csprefs.segment.output.CSF` - Option for writing the CSF. Default is `[0, 0, 0]`.

- f. `csprefs.segment.output.biascor` – Bias correction. A value of 1 means save bias corrected.
- g. `csprefs.segment.output.cleanup` - Clean up any partitions. Options are 0, 1 and 2. Each option is explained below.
  - a. 0 - Dont do cleanup
  - b. 1 - Light Clean
  - c. 2 - Thorough Clean

## **Regular Expressions**

Regular expressions are used to do pattern matching. Some of the standard regular expressions are as follows:

- a. `'\<\w+\>'` – String containing alphabets, numerals or underscore characters like `'ab_12_dd'`, `'sub01_vis'`, `'sub1'`, etc.
- b. `'\<d\>'` – String containing only numerals like `'1'`, `'2'`, etc.
- c. `'\<d{8}_d{6}_d{8}\>'` – String containing 8 decimals followed by underscore character followed by 6 decimals followed by underscore character followed by 8 decimals like `'20040331_111557_04030569'`, `'20040604_104632_04030708'`, etc.